

# Converting the PMT Container Testing Raw Data to ROOT File Format

Email: [zhaor25@mail2.sysu.edu.cn](mailto:zhaor25@mail2.sysu.edu.cn)

School of Physics



中山大學  
SUN YAT-SEN UNIVERSITY



# Outline

- ① converting raw data to root file
- ② update of container test results
- ③ Summary

# motivation

- 1 The Raw data of PMT testing is significant for the evaluation of PMT performance.
- 2 **While,Currently, the raw data of container system is not well organized and it is not convinient for people to get a quickly access.**
- 3 It is useful to convert all the testing raw data to ROOT format.
  - decrease the file size
  - easy to analysis and manage.
  - shadow the hardware details.

# requirements

- ① store the raw waveform data(.1pe, 1pe, TTS).
- ② store the auxiliary testing information(container , mass, HV, DCR. etc).
- ③ easy to manage (create, modify and update) and analyze.
- ④ one can acquire almost all the data needed for analysis(of one PMT) from only one file rather than collecting the details from server.

# prliminary file structure and strategies

- each PMT have one root file named in "SN\_rawdata.root"
- In a specific root file, we have several waveform waveform trees and a auxiliary data tree
- if one PMT go through several tests in the container, all the data will be saved still in only one root file but with different name of wave trees<sup>1</sup>; their auxiliary information will be filled several times in the same tree.





# current states

finished:

- basic structure of root file and TTree
- example rawdata root file
- example cpp program to access the waveforms from the generated root file

still working on :

- refine the root file contents and structure
- writing the document for potential users

# updates from shanghai colaberation meeting



# summary

- the converting of raw data from binary to root format is almost done
- one can easily restore the test waveforms with no loss of information
- the file size<sup>2</sup> decrease about 20% after transform<sup>3</sup>.
- the update of container results

---


<sup>2</sup>the total additional disk space requirement for 20k is less than 2T, so this not a problem

<sup>3</sup>about 50MB for one PMT of one light intensity

# THANKS

# BACK-UP

# TTS of HAMAMATSU PMT



hit-time and trigger time